

Marietta Memorial Hospital Department of Gastroenterology Juan Tejada, MD| Ben Scheinfeld, MD| Basile Njei MD Shelley Day, NP| Rosa Nardo, PA| Ashley Moore, NP| Ashley Wiggers, NP 400 Matthew Street, Suite 201 Marietta, OH 45750 | (740) 568-4590 | Fax (740) 568-4592 799 Farson Street, Suite 212 Belpre, Oh 45714

Helicobacter pylori Breath Test

Patient Preparation:

1. Patient should be fasting for 1 hour.

2. Do not administer this test if this list of instructions is not followed, as test results may be inaccurate:

a. Patients should not take proton-pump inhibitors (eg, Prilosec, Prevacid, Aciphex, Protonix, and Nexium) or bismuth compounds (eg, Pepto-Bismol) for 2 weeks prior to testing.

b. Patients should not take antibiotics for 4 weeks prior to testing.

c. Histamine 2-receptor antagonists (H2RA), such as Pepcid, Tagamet, Axid, or Zantac may impact the sensitivity of the test and, if possible, should be discontinued for 2 weeks prior to testing.

3. Carafate (sucralfate) does not interfere with the test. Use of antacids does not affect the accuracy of this assay.

Supplies: H. Pylori Breath Kit

Collection Instructions:

1. Do not collect if patient is younger than 3 years of age.

- 2. Follow instructions included with kit.
- 3. Mixing the (13)C-Urea Tablet

a. Dissolve the Citrica and the (13)C-enriched urea tablet in 150 to 200 mL (5.1 to 6.8 oz.) of tap water in the provided drinking cup.



Marietta Memorial Hospital Department of Gastroenterology Juan Tejada, MD| Ben Scheinfeld, MD| Basile Njei MD Shelley Day, NP| Rosa Nardo, PA| Ashley Moore, NP| Ashley Wiggers, NP 400 Matthew Street, Suite 201 Marietta, OH 45750 | (740) 568-4590 | Fax (740) 568-4592 799 Farson Street, Suite 212 Belpre, Oh 45714

b. Close the lid firmly using both hands. Place fingers over lid and shake thoroughly for a few minutes, until the Citrica Powder and the urea tablet are completely dissolved.

Note: Tiny particles may remain visible after thorough mixing. However, if more substantial particulate matter is still present after five minutes of mixing, discard the solution and repeat the procedure with a new kit.